

# Pick Your Poison. Isolation of Paclitaxel.

The diterpenoid, paclitaxel, which was identified as a potent chemotherapy agent for breast and ovarian cancer originates from the Pacific Yew tree. Though the life-saving compound was initially found in the bark, its presence was later confirmed in every part of the tree excluding the fruit.<sup>1</sup> Paclitaxel, is a member of the taxane family and has been associated with the highly toxic nature of the *Taxus* genus.<sup>2</sup> Paclitaxel was discovered to inhibit mitosis in tumor cells and other fast growing cells through the binding of microtubules.<sup>3</sup> This behavior led to its use as a tumor suppression agent. The isolation of paclitaxel from its major impurities is shown below with the use of Hamilton's PRP-1 (5 µm) HPLC column.

The isolation of paclitaxel was achieved by utilizing a sodium citrate buffer with a concentration of 20 mM at pH 6.0. The aqueous modifier promotes a stronger interaction between the analytes and the highly lipophilic resin. Complete baseline separation was observed using a 0.75 mL/minute flow rate at 35°C. While other HPLC methods employ a 250 x 4.6 mm HPLC column, we were able to reduce the column dimensions to 150 x 4.1 mm to facilitate a faster method while still maintaining good peak shape and resolution. Scale-up is easily accomplished due to Hamilton Company's proprietary stationary phase which can tolerate large sample loads without sacrificing peak shape, resolution or exhibiting carryover peaks.

- 1) Wilson, C., Hooser, S. *Toxicity of Yew Alkaloids, Veterinary Toxicology (Third Edition)*, Academic Press, 2018, 66, 947.
- 2) Muggia, F.; Kudlowitz, D. *Anti-Cancer Drugs*, 2014, 25 (5), 593.
- 3) Cortes, J.; Vidal, M. *Breast Cancer Res. Treat.* 2012, 133 (3), 821.


Author: Adam L. Moore, PhD.

## Column Information

<b>Packing Material</b>	PRP-1, 5 µm
<b>P/N</b>	79444

## Chromatographic Conditions

<b>Gradient</b>	Isocratic
<b>Temperature</b>	35°C
<b>Injection Volume</b>	5 µL
<b>Detection</b>	UV at 230 nm
<b>Dimensions</b>	150 x 4.1 mm
<b>Eluent A</b>	1:1 20 mM citrate pH 6:ACN
<b>Flow Rate</b>	0.75 mL/min.



©2021 Hamilton Company. All rights reserved.  
All other trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.  
Lit. No. L80119 — 10/2021