




96-Deep Well Storage Plates

Made in the U.S.A.




219009

- Manufactured from virgin polypropylene which has been pre-tested for low extractables
 - Alphanumeric grid-referencing
 - DNase / RNase free
 - Packed in sealed sleeves of 5 plates
- Rimmed version to stop cross contamination and enable a better seal
- Non-rimmed version to allow insertion of vials, or where the seal is not crucial
 - Cylindrical well with round bottom for optimal mixing and recovery
 - No inner edges to allow better collection of magnetic beads
 - Can be stored at -80°C
 - Plates are available non-sterile and sterile



219002



WITH RIM

ALL DIMENSIONS IN MM

| Catalog No. | Description | Case |
|-------------|--|------|
| 219002 | 96 Round-Well Storage Plate, 1mL, 42mm Height, Round Bottom, Raised Rimmed | 50 |
| 219037 | 96 Round-Well Storage Plate, 1mL, 42mm Height, Round Bottom, No Rim | 50 |
| 219009 | 96 Square-Well Storage Plate, 2mL, 44mm Height, Pyramid Bottom | 50 |