

Microlab 600 Basic Single Syringe Dispenser

Specification Sheet (Part/REF # ML610-DIS)



The Single Syringe Dispenser is ideal for dedicated dispensing applications requiring a single syringe drive. Improve accuracy, throughput, and consistency when compared to traditional laboratory equipment.

Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Program Memory	2 GB
Valve Setup	Single Syringe Dispenser
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder
Controller	Basic
Number of Syringe Drives	One
Probe	Concorde CT Hand Probe
Syringes Included	Includes one user selected syringe
Fluid Path Materials	FEP, CTFE, PTFE, Borosilicate Glass

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/microlab-600/ml610-dis>

Spec. Version B