

Microlab 600 Advanced Single Syringe Dispenser

Specification Sheet (Part/REF # ML620-DIS)



The Single Syringe Dispenser is ideal for dedicated dispensing applications requiring a single syringe drive. Improve accuracy, throughput, and consistency when compared to traditional laboratory equipment.

Product Specifications

| | |
|---------------------------------|---|
| Power Rating | 24 VDC, 2.5A |
| Power Requirements | 100-240V 1.5A max 50-60 Hz |
| Program Memory | 2 GB |
| Valve Setup | Single Syringe Dispenser |
| Syringe Drive Mechanism | Stepper motor driven high precision lead screw with encoder |
| Controller | Advanced |
| Number of Syringe Drives | One |
| Probe | Concorde CT Hand Probe |
| Syringes Included | Includes one user selected syringe |
| Fluid Path Materials | FEP, CTFE, PTFE, Borosilicate Glass |

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/microlab-600/ml620-dis>

Spec. Version B