

Microlab 600 Advanced Dual Syringe Dispenser

Specification Sheet (Part/REF # ML625-DIS)



The Dual Syringe Dispenser has two syringe drives for applications such as Animal Dosing that require the use of more than one liquid at the same time. The two liquids can be dosed from the same instrument saving setup and dosing time.

Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Program Memory	2 GB
Valve Setup	Dual Syringe Dispenser
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder
Controller	Advanced
Number of Syringe Drives	Two
Probe	Dual Push Button Hand Probe
Syringes Included	Includes two user selected syringes
Fluid Path Materials	FEP, CTFE, PTFE, Borosilicate Glass

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/microlab-600/ml625-dis>

Spec. Version B